

## Check-Out

- ☐ MSR geräumt (Liege steht neben den Dewar, Bildschirm und Spiegel in der Ausgangsposition)
- ☐ Dewar in 25° Stellung – *Positionierungsanzeige blinkt grün*
  
- ☐ MSR-Tür geschlossen
- ☐ Licht aus im MSR und Einstellung der Helligkeit bei 12:00 Uhr
- ☐ Stimulus-Cabinet geräumt (Steckdosenleiste aus)
- ☐ Lichtprojektor aus
- ☐ Licht aus im MEG-Raum
  
  
- ☐ Acquisition-PC – Sitzung geschlossen
- ☐ Acquisition-Bildschirm aus
- ☐ Überwachungsmonitor aus
- ☐ Kopf-Digitalisierungssystem aus
- ☐ Fenster geschlossen
- ☐ Licht aus im MEG-Vorraum
- ☐ Tür zum Flur abschließen

- Verlassen Sie die Räumlichkeiten wie Sie sie vorgefunden haben.

Bitte schließen Sie die Tür zum Gang ab.

## Check-Out

- ☐ MSR room cleared (store lounger, reset screen mirrors to default)
- ☐ Dewar position at 25 degrees – check lifting mechanism indicator flashes green
  
- ☐ MSR-Door closed
- ☐ Light OFF in the MSR and 12:00 o'clock position for brightness
- ☐ Stimulus-Cabinet cleared (power strip off)
- ☐ Projector - power-off
- ☐ Light off MEG-room
  
- ☐ Acquisition-PC - session closed
- ☐ Acquisition-Monitor off
- ☐ CCTV-Monitor off
- ☐ Head-registration system off
- ☐ Windows closed
- ☐ Light off MEG-antechamber
- ☐ Door closed

- Please leave MEG cleaner than you found it

Please lock the door to the aisle.